

# Hindustan Syringes & Medical Devices Ltd.

## INSTRUCTIONS FOR USE

### DISPOWAY® 3 Way Stop Cock

**Manufacturer:**

Hindustan Syringes & Medical Devices Ltd, 99/25 Ballabgarh, Faridabad, Haryana 121004 India.

**QUALITATIVE AND QUANTITATIVE COMPOSITION:**

3 Way stop Cock body & luer lock connector made of Poly carbonate, core made of Poly ethylene, luer slip & luer lock stopper are made of Polypropylene.

**INDICATIONS:**

- This medical device is used to connect one or more I.V. infusion set in order to aspirate fluids or to administer medicinal products in liquid form intracutaneously, subcutaneously, intramuscularly, intravenously.

**CONTRAINDICATIONS:**

- It must not be used in patients with known hypersensitivity to any of the components.
- Administration of highly viscous fluids.
- Massive blood transfusions.

**INSTRUCTIONS FOR USE:**

- Carefully examine the integrity of the product packaging and the expiry date;
- Open individual packaging at the marking and then remove this medical device from the packaging;
- Connect one or two I.V. infusion lines as required to the ports;
- Remove the luer slip cap from the male portion & hold it upright;
- Open clamp of I.V. infusion lines & displace air entirely by rotating the tap to open position;
- Rotate tap to close position;
- Connect it to venipuncture device & secure it with connector. Do not use excessive force while tightening.
- Rotate the tap to open position and adjust the flow of infusion set.

**STORAGE CONDITIONS:**

- The product should be kept in the original packaging, protected from humidity, excessive temperature and light, without corrosive / Acidic Fumes.
- Do not refrigerate. Do not freeze.

**WARNINGS:**















1. This product may be used only by a qualified medical doctor or nurse.
2. Before using the medical device please read carefully the Instructions for Use.
3. The product should be used in accordance with the Instructions for Use (see above).
4. For single use only. Reuse of Device may lead to cross contamination / Infection.
5. The product is sterilized by ethylene oxide. Do not resterilize. Do not recycle.
6. Before use, ensure that the packaging is not damaged. Irregular transportation, storage and handling may cause structural and/or functional damages to the product or its packaging. If there is any damage inflicted to the product do not use it. Use this medical device immediately after the opening of the individual packaging.
7. This product should be used immediately after opening it.
8. NON-TOXIC, STERILE, NON-PYROGENIC if the packaging is not opened or damaged.
9. It does not cause hemolysis.
10. Do not use after the expiry date.

**Shelf life:** 5 years.

**Date of the last revision of the text:**

Dec. 2023

**Symbols:**

	Manufacturing Date
	Expiration Date
	Reference Number
	Keep away from sunlight
	Non-pyrogenic
	Instruction for Use
	Do not re-sterilize
	Do not use if package is damaged
	Keep dry
	Do not re-use
	Medical Device
	Caution
	Sterilized using ethylene oxide
	Discard in a sharp container

**LIABILITY:**

Hindustan Syringes & Medical Devices Ltd. warrants that its products are free from defects in workmanship and material under proper handling and use. This warranty is void if a failure or a defect is due to improper handling or storage or factors related to the patient, his/her diagnosis, treatment, surgical procedures or any other matter beyond the control of the manufacturer Hindustan Syringes & Medical Devices Ltd. that has direct effects on the products of Hindustan Syringes & Medical Devices Ltd. and the results obtained from their use. The liability of Hindustan Syringes & Medical Devices Ltd. is limited up to the replacement of the product only. Hindustan Syringes & Medical Devices Ltd. disclaims any responsibility for possible consequences resulting from improper use or use in the procedures for which this medical device is not indicated.