




DISPO VAN®

HOW TO USE AN INSULIN SYRINGE



1.

Pull off the bottom cap and then the needle cover.
2.

Pull the plunger to draw air into the syringe, equal to the dose of insulin to be given.
3.

Penetrate needle into the center of the vial, push the plunger down to inject the air.
4.

Hold vial upside down, to draw the prescribed dose of the insulin while ensuring no air bubbles should be present.
5.

Disinfect chosen site, using an alcohol swab insert the needle straight (90 degree angle) into the skin as per the recommended sites.
6.

Push the plunger slowly and steadily until the insulin dose is completely injected.
7.

Pull out the needle from skin while maintaining a 90 degree angle.
8.

Dispose off the used syringes into the sharp container for safe disposal.

- If air bubbles appear in the syringe, push the plunger up, injecting insulin back into the vial redraw insulin to the desired number of units.
- To minimize the risk of needle stick injury, use safety box to dispose off used medical device.
- The use of a short needle without a skinfold is recommended for most patients when given by a healthcare worker or caregiver.
- Young patients of 6 year age and younger as well as very thin adults may require a pinch up, prior insertion.

SYRINGE CAPACITY	NEEDLE SIZE	LENGTH	INSULIN UNIT ADMINISTRATION	
			U-100	U-40
0.3 ml	29G X 1/2"	(0.33 X 13 mm)	30 units or less of insulin concentration U-100	-
	30G X 5/16"	(0.30 X 8 mm)		
	31G X 15/64"	(0.25 X 6 mm)		
	31G X 3/16"	(0.25 X 5 mm)		
0.5 ml	29G X 1/2"	(0.33 X 13 mm)	50 units or less of insulin concentration U-100	-
	30G X 5/16"	(0.30 X 8 mm)		
	31G X 15/64"	(0.25 X 6 mm)		
	31G X 3/16"	(0.25 X 5 mm)		
1 ml	29G X 1/2"	(0.33 X 13 mm)	100 units or less of insulin concentration U-100	40 units or less of insulin concentration U-40
	30G X 5/16"	(0.30 X 8 mm)		
	31G X 15/64"	(0.25 X 6 mm)		
	31G X 3/16"	(0.25 X 5 mm)		

CAUTION

- DO NOT ATTEMPT TO RESTERILIZE AND REUSE
- Before use check integrity of tamper proof heat staked seals where end caps and barrel meet. The fluid path of the syringe is guaranteed to be sterile if the syringe has not been opened or damaged.
- For patients who mix two types of insulin, do not change (a) the order of mixing that the physician has prescribed or (b) the model and brand of syringe and needle without first consulting your physician or pharmacist. Failure to heed this warning can result in dosage error.
- Before changing to a shorter needle consult your physician.
- RECOMMENDED FOR PERSONAL USE.

REPORTING OF SERIOUS ADVERSE EVENTS:

- In case of any serious adverse events, please report the same to E-mail: info@hmdhealthcare.com
- For product used in EU also report to the competent authority of the member state in which the user and/or patient is established.

	Caution
	Batch code
	Date of manufacture
	Use - by date
	Consult Instructions for use
	Do not re-use
	Do not use if package is damaged
	Do not re-sterilize
	Sterile fluid path Sterilized using ethylene oxide
	Medical Device
	Single sterile barrier system with protective packaging outside
	Non-pyrogenic
	Keep away from sunlight
	Keep dry
	Discard in a sharps container



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
174,178/25 BALLABGARH, FARIDABAD, HARYANA 121004 INDIA
M.L. MFG/MD/2020/000307
www.hmdhealthcare.com

EC REP

Advena Ltd
Tower Business Centre, 2nd Fl.,
Tower Street, Swatar, BKR 4013 Malta

