KIT KATH[†]

KIT KATH[†] I.V. CANNULA

Instructions for use :

EN	KIT KATH [†] I.V. CANNULA	NO	KIT KATH [†] L.V. ANNULA
DE	KIT KATH [†] I.V. KANÜLE	PL	KIT KATH [†] KANIULA DOŻYLNÁ
AR	القنية الوريدية (IV) KIT KATH [†]	PT	KIT KATH [†] CÂNULA I.V.
HR	KIT KATH [†] L.V. KANILA	RO	KIT KATH [†] CANULĂ INTRAVENOASĂ
CS	KIT KATH [†] L.V. KANYLA	RU	KIT KATH [†] ВНУТРИВЕННАЯ КАНЮЛЯ
DA	KIT KATH [†] I.V. KANYLE	SR	KIT KATH [†] Л.В . ЦАННУЛА
NL	KIT KATH [†] I.V. CANULE	ES	CÁNULA INTRAVENOSA KIT KATH [†]
FI	KIT KATH [†] I.V. KANYYLI	SV	KIT KATH [†] I.V. KANYL
FR	KIT KATH [†] CANULE INTRAVEINEUSE (I.V)	TR	KIT KATH [†] I.v. KANÜL
EL	KIT KATH [†] Λ.Β. ΚΑΝΟΥΛΑ	SQ	KIT KATH [†] I.V. CANNULA
IT	KIT KATH [†] I.V. CANNULA	SK	I. V. KANYLY KIT KATH [†]
LT	KIT KATH [†] INTRAVENINIS KATETERIS		

MANUFACTURING SITE-1
**HINDUSTAN SYRINGES
& MEDICAL DEVICES LTD.**
174, 178/25 BALLABGARH, FARIDABAD
HARYANA 121004 INDIA

MANUFACTURING SITE-2
HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
JHARSANTLY-JAJRU ROAD, SECTOR-59 PART-II, BALLABGARH,
FARIDABAD, HARYANA 121004 INDIA
Made in India

Advena Ltd.
Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

MAY 2025/REV.-2

EN

KIT KATH® I.V. CANNULA

INSTRUCTIONS FOR USE

DEVICE DESCRIPTION :

- IV cannula is a small, flexible tube designed to deliver IV medications and fluids to the patient through the peripheral vein situated in the arm or leg of the patient.

MATERIALS USED

- PP, LDPE, HDPE, PC, Stainless steel
- Master Batch : Orange, Medium Grey, White, Deep Green, Pink, Deep Blue, Yellow & Violet
- This product does not contain natural rubber latex. This product is DEHP free.

INTENDED USE/INDICATIONS:

- IV cannula is a conventional device used for intravenous / intravascular access for short term peripheral cannulation.

CONTRAINDICATIONS:

- Product should not be used in patients with known hypersensitivity to any of the materials used.

INTENDED USER:

- Use is restricted to a qualified doctor or a paramedic.

INTENDED PATIENT POPULATION:

- All patient population irrespective of age and gender.

INSTRUCTIONS FOR USE:

- Carefully select and aseptically prepare the site.
- Select suitable size of IV. Cannula & inspect visually to ascertain that package is intact.
- Remove cannula from sterile packing.
- Grip the cannula from injection port cap & projection provided on hub (for IV. Cannula having injection port).
- Grip the cannula so the needle and catheter tip are properly aligned (for IV. Cannula not having injection port).
- Remove the needle cover.
- Perform venipuncture & check for flash back of blood in flash back chamber.
- Advance the catheter into vein and simultaneously withdraw the needle.
- Withdraw the needle completely while pressing the vein just after the tip of the catheter into the vein & discard the needle in appropriate container. Secure the Cannula to the patient with the help of adhesive dressing.
- Connect to the IV. Infusion set line.
- Cover the puncture site with sterile dressing.
- Drugs can be injected with the help of syringe without needle through integrated injection port after removing the port cap. Close the injection port after use (for IV. Cannula having injection port).
- Perform routine monitoring & venipuncture site maintenance according to medical norms.

PRECAUTIONS :

- Open the sterile barrier pack only under aseptic conditions.
- Never try to reinsert the partially or completely withdrawn needle.
- Close the port cap and stoppers when the device is not in use.
- Never rotate the needle or catheter while insertion.
- Select and appropriate size of device before use.
- Insert the needle with the bevel facing up and angle between 10° to 30°.
- Observe the flash back of blood into needle hub and confirm the puncture of the vein.

CAUTION :

- The use of this product is restricted to a qualified doctor or a paramedic
- Read instructions before use.
- HMD DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.
- Store in clean, dry & insect free place.
- Avoid exposure to Acidic/Corrosive fumes to avoid corrosion/Rust.
- Do not expose to excessive heat or direct sunlight as it can affect packing/sealing integrity.
- Injection from port should be done at slow speed, avoid trying to inject rapidly and suddenly (for IV. Cannula having injection port).
- The product should be used immediately after opening the packaging.
- Visually inspect and carefully check the product and packaging before use. Improper transport, handling and storage may cause corrosion, structural and/or functional damage to device or packaging.
- Do not attempt to reinsert a partially or completely withdrawn needle, to avoid puncturing of catheter.
- Discard after use in sharps container.
- The product should not be reprocessed.
- For single use only. Do not clean or resterilize. Reuse of device may lead to cross contamination/ infection.
- The product is guaranteed non-toxic , sterile & non-pyrogenic if the package has not been opened or damaged.

ADVERSE EFFECTS:

Pain, discomfort, hives, bleeding, blanching, bruising, burning, allergic reaction from hypersensitivity to materials used in the device, local insertion site infection etc.

DISPOSAL:

After use, this product may be a potential biohazard. Handle and dispose off in accordance with accepted medical practice and applicable local, state, and federal law and regulation.

REPORTING OF SERIOUS ADVERSE EVENTS:

- In case of any serious adverse events, please report the same to E-mail : info@hmdhealthcare.com
- For product used in EU also report to the competent authority of the Member State in which the user and/or patient is established

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