
**DISPO VAN[®]**

HOW TO USE AN INSULIN SYRINGE

Instructions for use :

S.No.	LANGUAGE	TRANSLATION
1	EN (ENGLISH)	Single Use Insulin Syringe
2	HR (CROATIAN)	Insulinska štrcaljka za jednokratnu upotrebu
3	EL (GREEK)	Σύριγγα ινσουλίνης μιας χρήσης
4	ET (ESTONIAN)	Ühekordne insuliini süstal
5	FR (FRENCH)	Seringue à insuline à usage unique
6	DE (GERMAN)	Insulinspritze zum einmaligen Gebrauch
7	IT (ITALIAN)	Siringa da insulina monouso
8	LT (LITHUANIAN)	Vienkartinis insulino švirkštas
9	NL (DUTCH)	Insulinespuit voor eenmalig gebruik
10	PL (POLISH)	Jednorazowa strzykawka insulinowa
11	SK (SLOVAK)	Insulinová striekačka na jednorazové použitie
12	ES (SPANISH)	Jeringa de insulina de un solo uso
13	TR (TURKISH)	Tek kullanımlık İnsülin Şiringası
14	DA (DANISH)	Insulinsprøjte til engangsbrug

HINDUSTAN SYRINGS & MEDICAL DEVICES LTD.
17A, 17B/25 BALLABGARH, FARIDABAD, HARYANA 121004 INDIA
M.L. 9961602001000007
www.hmdhealthcare.com

Made in India

Advana Ltd
Tower Business Centre,
2 - Fl., Tower Street, Swatara,
BKR 4013 Malta

0297

AUG 2009/REV -0

EN

- If air bubbles appear in the syringe, push the plunger up, expelling insulin back into the vial (return insulin to the desired number of units).
- To minimise the risk of needle stick injury, use safety box to dispose of used medical devices.
- The use of a short needle without a blank is recommended for most patients when given by a healthcare worker or caregiver.
- Young children of 6 year age and younger as well as very thin adults may require a pinch up, prior injection.

Variant Description

SYRINGE CAPACITY	NEEDLE SIZE	LENGTH	INSULIN UNIT ADMINISTRATION	
			U-100	U-40
0.3 ml	29G X 1/2"	(0.33 X 13 mm)	30 units or less of insulin concentration U-100	
	30G X 5/16"	(0.30 X 8 mm)		
	31G X 15/64"	(0.25 X 6 mm)		
	31G X 3/16"	(0.25 X 5 mm)		
0.5 ml	29G X 1/2"	(0.33 X 13 mm)	50 units or less of insulin concentration U-100	
	30G X 5/16"	(0.30 X 8 mm)		
	31G X 15/64"	(0.25 X 6 mm)		
	31G X 3/16"	(0.25 X 5 mm)		
1 ml	29G X 1/2"	(0.33 X 13 mm)	100 units or less of insulin concentration U-100	40 units or less of insulin concentration U-40
	30G X 5/16"	(0.30 X 8 mm)		
	31G X 15/64"	(0.25 X 6 mm)		
	31G X 3/16"	(0.25 X 5 mm)		

CAUTION

- DO NOT ATTEMPT TO RESTERILIZE AND REUSE.
- Before use check integrity of tamper proof heat sealed seals where end caps and barrel meet. The fluid path of the syringe is guaranteed to be sterile if the syringe has not been opened or damaged.
- For patients who mix two types of insulin, do not change (a) the order of mixing that the physician has prescribed or (b) the model and brand of syringe and needle without first consulting your physician or pharmacist. Failure to heed this warning can result in dosage error.
- Before changing to a shorter needle consult your physician.
- RECOMMENDED FOR PERSONAL USE.

STORAGE CONDITIONS:

- Store in clean, dry and insect free place.
- Avoid exposure to excessive heat.
- Avoid exposure to acidic / corrosive fumes to avoid corrosion / rust.

REPORTING OF SERIOUS ADVERSE EVENTS:

- In case of any serious adverse events, please report the same to E-mail: info@finchhealthcare.com
- For product used in EU also report to the competent authority of the member state in which the user and/or patient is established.

SHELF LIFE: 5 years from the date of manufacturing.

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Symbol glossary

(symbol harmonizing with ISO 15223-1 : 2021)

	Caution
	Batch code
	Date of manufacture
	Use - by date
	Consult instructions for use
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Sterile fluid path sterilized using ethylene oxide
	Medical Device
	Single sterile barrier system with protective packaging outside
	Non-pyrogenic
	Keep away from sunlight
	Keep dry
	Discard in a sharps container
	Authorized Representative in the European Community
	Manufacturer
	Marking for meeting the requirements of European legislation The number displayed on the bottom is the notified body number
	Catalogue Number
	Unique Device Identifier
	Move vertically
	Fragile, handle with care

NOTE : Symbols used on the label and in this IFU apply only to this product and are based on current applicable standards at the time of release.