

Safe Injection



KOJAK SELINGE Auto Disable syringe for patient safety in the curative sector, a syringe that breaks automatically if you try to re-use it.

COMPONENTS:

BARRELS & PLUNGERS : are made of medical grade polypropylene that is non-toxic & compatible with the medication for short term contact drug delivery application.

GASKETS : Made as per standard IS:10258, chemically inert and has short term compatibility with the medication for short term contact drug delivery application.

LOCK : Presence of ring & lock mechanism which automatically breaks the plunger, if reuse is attempted.

NEEDLE : The tip of the needle has a thin wall and 3 facet lancet point. Needle is well siliconized to reduce friction (for smooth penetration & withdrawal).

PRECAUTIONS : Before use, check the integrity of the sealing of the pack. Do not push the plunger forward before withdrawing medicament from the vial otherwise plunger will get locked.

KEY BENEFITS OF KOJAK SELINGE :

- Significantly helps reduce infections due to unsafe injections.
- Reduces health-care costs at a macro-level.
- Reduce average stay of patient in hospital.
- Usage of fixed needle syringes assures that neither syringes nor the needle can be used again.

FOR CURATIVE

PRODUCT RANGE			
SIZE	NEEDLE SIZE		WHO PQS No.
2 ml	26G x 1/2"	24G x 1"	E013/001
3ml	*0.45 x 13	0.55 x 25	E013/018
5 ml	23G x 1"	22G x 1"	
	*0.60 x 25	0.70 x 25	
10 ml	22G x 1 1/4"	21G x 1 1/2"	E013/020
20 ml	*0.70 x 32	0.80 x 38	E013/039

Conform to ISO7886-4

* All units are in mm.

Auto Disable Syringe with/without - Fixed/Detachable Needle.