

# Safe Injection



**KOJAK SELINGE** Auto Disable syringe for patient safety in the immunization sector, a syringe that breaks automatically if you try to re-use it.

## COMPONENTS:

**BARRELS & PLUNGERS** : are made of medical grade polypropylene that is non-toxic & compatible with the medication for short term contact drug delivery application.

**GASKETS** : Made as per standard IS:10258, chemically inert and has short term compatibility with the medication for short term contact drug delivery application.

**LOCK** : Presence of ring & lock mechanism which automatically breaks the plunger, if reuse is attempted.

**NEEDLE** : The tip of the needle has a thin wall and 3 facet lancet point. Needle is well siliconized to reduce friction (for smooth penetration & withdrawal).

**PRECAUTIONS** : Before use, check the integrity of the sealing of the pack. Do not push the plunger forward before withdrawing medicament from the vial otherwise plunger will get locked.

## KEY BENEFITS OF KOJAK SELINGE :

- Significantly helps reduce infections due to unsafe injections.
- Reduces health-care costs at a macro-level.
- Reduce average stay of patient in hospital.
- Usage of fixed needle syringes assures that neither syringes nor the needle can be used again.

## FOR IMMUNIZATION

PRODUCT RANGE			
SIZE	NEEDLE SIZE		WHO PQS No.
0.1ml	26G x 3/8 0.45 x 10	27G x 3/8 0.40 x 10	E008/045
0.5 ml	23G x 1" 0.60 x 25		E008/039
1ml	26G x 1/2" * 0.45 x 13 23G x 1" *0.60 x 25	24G x 1" 0.55 x 25 22G x 1" 0.70 x 25	E013/056

Conform to ISO7886-3, IS10258-3

\* All units are in mm.

Auto Disable Syringe with Fixed Needle.