

PERFEKTUM™



PERFEKTUM™

Sterile Single Use Syringes

“Ribbon” Packed PERFEKTUM, single use syringes without needle & needle mounted on top, are synonymous with quality

COMPONENTS :

BARRELS & PLUNGERS : made of non-toxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application.

The nozzle of the barrel is of two types - luer mount and luer lock. Both have a 6% luer taper as per ISO:594.

GASKETS : Chemically inert and compatible with the medications for short term contact drug delivery application.

PACKING : PERFEKTUM, single use syringes are individually ribbon packed in a double laminated plastic film. Ribbon packaging is superior to blister packaging due to higher puncture and water resistance, enabling extended shelf-life under adverse conditions.

PRODUCT RANGE			
Capacity	Type of Nozzle		Standard Needle Dimensions
2 ml.**	Luer mount/ central nozzle	with	24G×1", 23G×1", 23G×1½", 22G×1" *0.55×25 0.60×25 0.60×38 0.70×25 22G×1¼", 22G×1½", 21G×1", 21G×1½" *0.70×32 0.70×38 0.80×25 0.80×38
3 ml.	Luer mount/ central nozzle	with	Specifications same as 2 ml
5 ml.**	Luer mount/ central nozzle	with	Specifications same as 2 ml
10 ml.**	Luer mount with eccentric nozzle		22G×1¼", 22G×1½", 21G×1½", 20G×1½" *0.70×32 0.70×38 0.80×38 0.90×38 Same as above
20 ml.	Luer mount with eccentric nozzle		Specifications same as 10 ml
50 ml.	Luer mount with eccentric nozzle		Without needle



**also available in Blister pack.

*All units are in mm

• Conform to : ISO:7886-1/IS:10258 standards.

▲ Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.



www.hmdhealthcare.com