

KOJAK

SAFETY BOX

For Collection of Sharps-used syringes/needles/blades, 5L/10L



As per WHO standard WHO / PQS / E10 / SB01.1

For Resistance to piercing, Tamper proof, Water resistance, Temperature resistance.
Ease of formation on reading pictorial instructions on the Safety Box.

SUPERIOR THICKNESS AND BOARD DENSITY.

SPECIFICATION:

Volume : 5ltr/10ltr.

CAPACITY:

S. No.	Size of Syringe	With Needle No. of Syringes without Agitation		*Without Needle No. of Syringes without Agitation		With Needle No. of Syringes with Agitation		*Without Needle No. of Syringes with Agitation	
		5 ltr	10 ltr	5 ltr	10 ltr	5ltr	10ltr	5ltr	10ltr
1.	0.5ml/2ml	100	200	275	630	150	275	295	650
2.	5 ml	65	120	140	320	120	190	170	360

* Defanged/Cut Needle (5 ltr - suitable for single use syringes 1 ml to 10 ml; 10 ltr suitable for single use syringes 1 ml to 20 ml)

Material	- Yellow Duplex Kraft Board
Thickness Of Walls	- Over 1.0mm
Grammage	- Over 850gsm

Technical Dimensions

Weight fully assembled	275gms (min) (5 ltr)	450gms (min) (10 ltr)
Diameter of syringe insert hole	38 mm (5 ltr)	42 mm (10 ltr)
Cartons	5 Ltr. (25 pcs)	10 Ltr. (20 pcs)
Shipping Weight	5 Ltr. - 7.5 kg(min.)	10 Ltr. - 10 kg(min.)

External size

5 Ltr- Before Assembly	(580 to 590) x (280 to 290) x (3 to 5) mm
After Assembly	285 x 160 x 125 \pm 5 mm
10 Ltr Before Assembly	(660 to 670) x (400 to 410) x (3 to 5) mm
After Assembly	285 x 255 x 160 \pm 5 mm

FEATURES/BENEFITS:

- Easy Assembly without training
- Direction of use on the box in multilingual text as well as pictorial directions and sequenced code on flaps.
- Handle in center avoids spillage/tilted Box during transportation.
- For Disposal - As per country specific waste disposal national policy guidelines.
- Seal edges with Adhesive Surgical Tape in case Box is intended to be used to collect Surgical Blades.



www.hmdhealthcare.com