

UNOLOK[®]

INSULIN SYRINGE

0.3 ml/0.5ml/1ml

The new measure of comfort for diabetic patients



- Choice of 3 syringe sizes with clear markings to ensure precise dosage.
- Lubricated barrels for smooth plunger movement.
- Utmost injection comfort, a proven quality of HMD syringes.

**New 31 G
Insulin
Syringe**



PRODUCT SPECIFICATIONS

PRODUCT

UNOLOK[®] INSULIN SYRINGE

PRODUCT RANGE

SPECIFICATIONS

- 0.3ml/0.5ml/1ml syringe with 31G x 15/64 , 30G X 5/16" & 29G X ½" needle.
- Shorter 31G needle helps in preventing the much thinner needle from bending.
- Twice accuracy of dosage in 0.3 / 0.5 ml syringe due to scale increments of 1 unit as compared to 2 units in 1 ml U 100 syringe.
- 1 ml is available in both U 40 and U 100. 0.3 ml and 0.5 ml are available in U 100.
- Available in economy multipacks of 10 pieces in a double laminated plastic film or single unit in a ribbon pack.
- Needle cover and bottom cap is tamper proof heat sealed with the barrel which ensures that the syringes interior is sterile. (In case of multipacks)
- Siliconised fine needle for minimal pain.
- Negligible dead space (<0.01 ml).
- Needle is integrated with the barrel
- Gasket is chemically insert and compatible with the medication for short term contact drug delivery application.
- Manufactured as per ISO:8537

PRECAUTIONS

- Before use check the integrity of tamper proof heat-stacked sealing where the end-caps and barrel meet (multi pack). Twist off the bottom cap and the needle cover to expose the plunger and the needle respectively.
- Pull straight off to avoid damaging needle point.

WARNING

- Do not share syringes to avoid cross contamination.
- Do not attempt to re-sterilise or reuse.

NOTE : Insist on **UNOLOK**[®] only. All single use syringes are not **UNOLOK**[®]

UNOLOK INSULIN SYRINGES

DESCRIPTION	SIZE	NO.OF PCS. PER CASE	CASE DIMENSION	WEIGHT PER CASE (1ml)
MULTI PACK	0.25 X 8mm / 0.25 X 6mm (31G X 5/16") / (31G X 15/64")	2000	770 X 390 X 195mm	10kg
UNIT PACK	0.25 X 8mm / 0.25 X 6mm (31G X 5/16") / (31G X 15/64")	2000	735 X 375 X 275mm	10kg
MULTI PACK	0.30 X 8mm (30G X 5/16")	2000	770 X 390 X 195mm	10kg
UNIT PACK	0.30 X 8mm (30G X 5/16")	2000	735 X 375 X 275mm	10kg
MULTI PACK	0.33 X 13mm (29G X ½")	2000	770 X 390 X 195mm	10kg
UNIT PACK	0.33 X 13mm (29G X ½")	2000	735 X 375 X 275mm	10kg

*For 0.3 ml/0.5ml weight is about 8 kg



www.hmdhealthcare.com