



UNOLOK[®]

Sterile Single Use Syringes

“Ribbon” Packed UNOLOK, single use syringes without needle & needle mounted on top, are synonymous with quality

COMPONENTS :

BARRELS & PLUNGERS : made of non-toxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application.


The nozzle of the barrel is of two types - luer mount and luer lock. Both have a 6% luer taper as per ISO:594.

GASKETS : Chemically inert and compatible with the medications for short term contact drug delivery application.

LUER LOCK nozzle present in UNOLOK syringes ensures a secure needle fitment.

PACKING : UNOLOK single use syringes are individually ribbon packed in a double laminated plastic film. Ribbon packaging is superior to blister packaging due to higher puncture and water resistance, enabling extended shelf-life under adverse conditions.

PRODUCT RANGE		
Capacity	Type of Nozzle	Standard Needle Dimensions
1 ml.* (T.B.)	Luer mount with central nozzle	26G × 1/2" *0.45 × 13
2 ml.**	Luer mount/luer lock with central nozzle	24G × 1", 23G × 1" 23G × 1 1/2", 22G × 1" *0.55 × 25 0.60 × 25 0.60 × 38 0.70 × 25 22G × 1 1/4", 22G × 1 1/2", 21G × 1", 21G × 1 1/2" *0.70 × 32 0.70 × 38 0.80 × 25 0.80 × 38
3 ml.	Luer mount/luer lock with central nozzle	Specifications same as 2 ml
5 ml.**	Luer mount/luer lock with central nozzle	Specifications same as 2 ml
10 ml.**	Luer mount with eccentric nozzle	22G × 1 1/4", 22G × 1 1/2", 21G × 1 1/2", 20G × 1 1/2" *0.70 × 32 0.70 × 38 0.80 × 38 0.90 × 38
	Luer lock with central nozzle	Same as above
20 ml.	Luer mount with eccentric nozzle Luer lock with central nozzle	Specifications same as 10 ml
50 ml.	Luer mount with eccentric nozzle Luer lock with central nozzle	Without needle



**also available in Blister pack.

*All units are in mm

• Conform to : ISO:7886-1/IS:10258 standards.

▲ Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

