

# KOJAK<sup>®</sup>

## Disposable Sharps Container

For safe disposal of used syringes, needles & sharps



# DESCRIPTION

1. Suitable for handling used Syringes, Needles & Sharps.
2. Anti-puncture Durable plastic shell, made of polypropylene.
3. Comes with an easy to use rotatable lid
4. Available in 2L and 5L.
5. Designed to facilitate safe disposal to reduce NSI's (Needle Stick Injuries)



Front



Top

VARIANT	HEIGHT	WIDTH	WEIGHT
2L	15 cm	13.5 cm	77 g
5L	20.5 cm	19.6 cm	171 g

## Caution

1. Appropriate care should be used when filling, transporting and disposing in order to avoid the possibility of injury.
2. Dispose as per waste disposal national policy guidelines.

## Warning

Do not overfull more than used capacity marking

