

VAKU-8TM

Blood Collection Needle



MDSAP
MEDICAL DEVICE SINGLE AUDIT PROGRAM



HINDUSTAN SPRINGERS
& MEDICAL DEVICES LTD.
ADDRESS ISO 13485 : 2016

PRODUCT SPECIFICATIONS

Product - "VAKU-8" Brand Multi sample Blood Collection Needle.

COMPONENTS

3.1 TOP NEEDLE COVER - Top cover are made of medical grade PP. It gives protection to needle from damage.

3.2 LOWER NEEDLE COVER - Lower covers are made of medical grade HDPE.

3.3 NEEDLE TUBE - Needle tube is made from cold roll stainless steel. The needle tube is thin walled which allows a larger internal diameter resulting in improved flow rate.

3.4 NEEDLE SHEATH - Flexible rubber tube covers the backside of needle tube. This act as end blood stopper when needle is withdrawn from sample tube.

3.5 SILICONISED NEEDLE - To minimize trauma and facilitate comfortable insertion.

3.6 SHORT BEVEL - To minimize risk of cross perforation of vein.

PACKAGING

UNIT (BLOOD COLLECTION NEEDLE)

Each unit is sealed by adhesive tamper evident label.

i. Intermediate Carton - 10 shelf carton are packed in a three ply corrugated carton.

ii. Shipping Carton - 5 intermediate cartons are packed in a five ply corrugated carton.

Shipping Carton Size - 435 X 395 X 270 mm.

PRECAUTIONS

Before use check for sterility assurance and tamper evident sealing.

METHOD OF STERILIZATION

The product is sterilized by the ethylene oxide gas (EO₂).

PRODUCT RANGE

Needle Gauge	Colour Code	Needle Length (mm)	Needle Length (Inches)	Needle Diameter (mm)
21G	Green	38	1 1/2"	0.8
22G	Black	38	1 1/2"	0.7
23G	Blue	38	1 1/2"	0.6
21G	Green	25	1"	0.8
22G	Black	25	1"	0.7
23G	Blue	25	1"	0.6

WARNING

Do not use if the pack is previously opened or damaged. Ensure that blood collection tube is fully pressed over the blood collection needle. Do not attempt to re-sterilize and reuse.

Blood Collection Needle should not be recycled and should be preferably disposed off in sharps containers.

SHELF LIFE

The shelf life of product is 3 years from the date of manufacturing.

REFERENCE

Company's specification.



Scan to
Visit Product Page



Scan to
Visit Product Video