



**Cathy<sup>®</sup>** 

**Safety I.V. Cannula**  
with **SIPCLIP<sup>™</sup>**

*Taking care of those who take care of us*



*As part of your job and your passion, you care for hundreds of patients every day. Doctor, nurses and medical staff are depended upon in times of need. As a health professional you can now depend on someone as well.*

## HEALTHCARE PROFESSIONALS ARE AT HIGH RISK

As per WHO, among the 35 million health professionals worldwide, about 3 million receive percutaneous exposures to blood borne pathogens each year, 90% of these infections occur in developing countries.

## WHENEVER YOU ARE ON DUTY, CATHY IS ALSO ON DUTY

A protection device **SIPCLIP™** encapsulates the tip of the used needle when introducer needle is extracted from the catheter; Post Cannulation. SIPCLIP prevents accidental needle stick injury provided that general effective and safe working procedures and precautions are maintained during its use and disposal.

## AUTO PROTECTION

The hassel-free and auto protection technology of CATHY I.V Cannula is designed to provide “Involuntary Activated Safety” - an intuitive mechanism that requires no additional training of the health staff and strives to make I.V insertion and introducer needle extraction completely safe.

# COMPONENTS

### NEEDLE

Siliconised electrofinished, back cut bevelled stainless steel needle  
Instant Flash Back Notch on 24G and 26G Cannula

**SIPCLIP™**  
Sharp Injury  
Prevention Clip.



### LUER LOCK CAP

Closure to wing housing permits patient mobility when infusion temporarily discontinued.

### BLOOD STOPPER

The stopper has fine air vent to expedite blood flash back for enabling correct cannulation.

### NEEDLE HUB

Allows blood flash back visualisation.  
Thumb rest for proper grip during insertion.

### CATHETER

Available in VOLEX™ polyuerthane radio detectable.

### INJECTION PORT CAP

It facilitates extra medication and also prevents its back flow.

### BLISTER PACKAGE

Cannula is packed in rigid blister and sealed in peel open type medical grade

## PRODUCT RANGE

Colour Code	Gauge	Catheter Ext. Dia x Length (mm)	Water flow-rate (ml/min)
Grey	16G	1.70 x 45	200
White	17G	1.50 x 45	140
Green	18G	1.20 x 38	105
		1.20 x 45	100
Pink	20G	1.00 x 32	64
Blue	22G	0.80 x 25	38



Scan to  
Visit Product Video



Scan to  
Visit Product Page