



DISPO VAN[®]

INSULIN SYRINGE

TO ADMINISTER MEDICATION, COMFORTABLY

PRODUCT SPECIFICATIONS

Product Range

- 0.3ml/0.5ml/1 ml syringe with 31 G x 5/64, 30G X 5/16" & 29G X ½" needle.
- Shorter 31 G needle helps in dispersion of medicine in subcutaneous layer of skin & preventing needle from bending.
- Twice accuracy of dosage in 0.3 / 0.5 ml syringe due to scale increments of 1 unit as compared to 2 units in 1 ml U 100 syringe.
- Available in economical multipack of 10 pieces in a double laminated plastic film or as single unit in both ribbon pack and blister pack.
- Needle cover and bottom cap is tamper proof, heat sealed with the barrel ensuring sterility in the syringes . (In case of multipacks)
- Siliconised fine needle for comfortable insertion.
- Negligible dead space saves medication with no wastage (<0.01 ml).
- (Manufactured as per ISO:8537)



Product Range

Needle Size	Colour code		Description
29G x 1/2" 0.33 x 13mm	U-40	RED	Multipack
	U-100	ORANGE	
29G x 1/2" 0.33 x 13mm	U-40	RED	Unipack
	U-100	ORANGE	
30G x 5/16" 0.30 x 8mm	U-40	RED	Multipack
	U-100	ORANGE	
30G x 5/16" 0.30 x 8mm	U-40	RED	Unipack
	U-100	ORANGE	
*31G x 15/64" 0.25 x 6mm	U-40	RED	Multipack
	U-100	ORANGE	
*31G x 15/64" 0.25 x 6mm	U-40	RED	Unipack
	U-100	ORANGE	

Warning

- Do not share syringes to avoid cross contamination.
- Do not attempt to re-sterilise or reuse.

Precautions

- Before use check the integrity of tamper proof heat-stacked sealing where the end-caps and barrel meet (multi pack). Twist off the bottom cap and the needle cover to expose the plunger and the needle respectively.
- Pull straight off to avoid damaging needle point.
- Manufactured as per ISO:8537



Scan to
Visit Product Page