



KOJAK[®] SELINGE

THE KOJAK AUTO-DISABLE HYPODERMIC SYRINGE

To cater to your need for a safe injection



COMPONENTS

BARRELS & PLUNGERS

They are made of medical grade polypropylene that is non-toxic & compatible with the medication for short term contact drug delivery application.

LOCK:

Presence of ring and lock mechanism which automatically breaks the plunger, if reuse is attempted.

PRECAUTIONS :

Before use, check the integrity of the sealing of the pack. Do not push the plunger forward before withdrawing medicament from the vial otherwise plunger will get locked.

GASKETS

They are made from natural rubber, are chemically inert & compatible with all kinds of medications for short term contact drug delivery application.

NEEDLE

The tip of the needle has a thin wall and 3 facet lancet point. Needle is well siliconized to reduce friction (for smooth penetration & withdrawal).

FEATURES

Thin-walled needles with 3 facet lancet point for smooth and minimal pain penetration and withdrawal

Chemically inert gaskets suitable for short-term contact delivery application

Polypropylene syringes with markings for proper measurements

Lock and break mechanism to automatically break the plunger on its first use

Fixed needle mechm for ensuring needle not reused

Available for both immunization and curative segments



KEY BENEFITS OF KOJAK SELINGE:

Mitigates the risk of cross-contamination and blood borne illnesses

The fixed needle ensures that the needles can't be used again.

Minimises healthcare costs at macro-level

Available in detachable and fixed needle variants

It can also be used for mixing drugs and medications

Available for both curative and immunisation use

Capacity ranges from 0.1 ml to 20ml

FOR IMMUNISATION PRODUCT RANGE

SIZE	NEEDLE		WHO PQS NO.
0.1 ml	26G x 3/8 0.45x10	26G x 3/8 0.45x10	E008/045
0.5 ml	23G x 1" 0.60 x 25		E008/039
1 ml	26G x 1/2" *0.45x13 23G x 1" *0.60x25	24G x 1" 0.55x25 22G x 1" 0.70x25	E008/056

Conform to ISO7886-4

* All units are in mm.

Auto Disable Syringe with/without - Fixed/
Detachable Needle.

FOR CURATIVE PRODUCT RANGE

SIZE	NEEDLE	WHO PQS NO.
2 ml	26G x 1" *0.45x13	E013/001
3 ml	24G x 1" 0.55x25	
0.5 ml	23G x 1" 22G x 1" 0.60 x 25 0.70 x 25	E013/018
10 ml	22G x 1 1/4" *0.70x32	E008/056
20 ml	21G x 1 1/2" *0.80x38	E013/039

Conform to ISO7886-3, IS I 0258-3

* All units are in mm.

Auto Disable Syringe with Fixed Needle.



Scan to
Visit Product Video



Scan to
Visit Product Page