



®

GLASS VAN

Surgical Blades

Improved Sharpness for added precision



GLASS VAN[®] Surgical Blades

- Choice of Surgeons, Dentists, Chiropodists
- Made from premium quality Swiss High Carbon Steel
- Manufactured as per ISO : 7740
- CE Marked under Medical Device Directive (93/42/EEC)
- Individually packed in an easy peel aluminium pack
- Additionally protected for corrosion by VCI coated paper
- Sterilized by Gamma Radiation of 2.5 Mrad to ensure sterility assurance level of 10⁻⁶

Product range

Category	Commerical Designation No.
For small-size fittings (Handle No. 3)	10, 11, 12, 15
For large-size fittings (Handle No. 4)	20, 21, 22, 23, 24, 25



Packaging dimensions

Size of Shipping Carton	Primary Packaging (pcs)	Shipping Carton (pcs)	Gross Weight (Kg.)
540 x 180 x 180	100	5000	9.00
565 x 170 x 170	50	5000	9.00



Scan to
Visit Product