


**techno**  
**cut**

# Scalpel & Surgical Blades

*The edge of Perfection*



# TECHNOCUT Surgical Blades

- ◆ Ideal Choice for Surgeons, Dentists, Chiropodists
- ◆ Made from premium quality Swiss Carbon Steel & Stainless Steel
- ◆ Manufactured as per ISO : 7740
- ◆ CE Marked under Medical Device Directive (93/42/EEC)
- ◆ US FDA Regulated
- ◆ Individually packed in an aluminium peel pack which is protected by VCI coated paper Sterility with Gamma Radiation of 2.5 Mrad to ensure sterility assurance level of 10<sup>-6</sup> whilst preventing corrosion

## Product Range

Category	Commercial Designation No.
* For small-size fittings (Handle No.3)	10, 10A, 11, 12, 12D, 13, 15, 15C
* For large-size fittings (Handle No.4)	18, 19, 20, 21, 22, 23, 24, 25
Miniature Blades	6100, 6200, 6400, 6700, Minimal Incision Blade No. 67
Podiatry Blades	312, 313, 314, 316
Stitch Cutters	Stitch Cutter



## Packing Dimensions

Size of Shipping Carton (mm)	Primary Packaging (pcs)	Shipping Carton (pcs)	Gross Weight (Kg.)	Volumetric Weight (Kg.)
*540 x 180 x 180	100	5000	10.00	3.00

# TECHNOCUT

## Surgical Scalpel with Plastic Handle

- ◆ Disposable Scalpel with high quality Swiss Stainless Steel Blades
- ◆ Ergonomically designed handle for better feel & precision
- ◆ Ethylene Oxide (EO) Gas sterilized
- ◆ Protected edge with easy to remove plastic sheath
- ◆ Easy Peel Blister Pack in Medical Grade Paper complying with EN 868

### Packing Dimensions

Size of Shipping Carton (mm)	Primary Packaging (pcs)	Shipping Carton (pcs)	Gross Weight (Kg.)	Volumetric Weight (Kg.)
545 x 225 x 470	10	1000	12.20	9.50



Scan to  
Visit Product Page