



DISPO VAN[®]

STERILE SINGLE USE NEEDLES



STERILE SINGLE USE NEEDLE

Manufactured using premium-quality cold-rolled stainless steel strips with the latest Japanese technology, Dispovan needles are ground on high-precision machines: assembled and packed on totally automated machines. Vision cameras ensure 100% automated quality control, enabling superior sharpness to minimize pain.

COMPONENTS

NEEDLE HUBS

Made of non-toxic, medical-grade Polypropylene compatible with the injectable medication.

CANNULA

BURR FREE - 3 FACET LANCET POINT are made from SUS 304 stainless steel

PACKING

The needles are packed in tamper-evident blister packs having easily peelable, medical-grade paper which ensure reliability of seal and product sterility.

PRECAUTIONS

Before use, check the integrity of the sealing of the pack. Store in clean Dry & insect free place. Avoid exposure to excessive heat. Avoid exposure to acidic/corrosive fumes, to avoid corrosion/rust.

Avoid piercing the needle through the wall in hard plastic container, it may result in needle blockage.

COLOUR CODE: As per ISO:6009.



SINGLE USE HYPODERMIC NEEDLE

COLOUR CODE AND SPECIFICATION

SIZE		COLOUR		LENGTH						
OD (mm)	Guage	Colour & Code		3/8"	1/2"	5/8"	3/4"	1"	1 1/4 "	1 1/2 "
0.40	27G		Grey							
0.45	26G		Brown							
0.50	25G		Orange							
0.55	24G		Medium purple							
0.60	23G		Blue							
0.70	22G		Black							
0.80	21G		Green							
0.90	20G		Yellow							
1.10	19G		Creamish							
1.20	18G		Pink							
1.60	16G		White							

Conform to : ISO: 7864/IS: I 0654 standards.

Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.



Scan to
Visit Product Page