

UNOLOK[®]

STERILE SINGLE USE NEEDLES



UNOLOK[®]

STERILE SINGLE USE NEEDLES

Manufactured using premium-quality cold-rolled stainless steel strips with the latest Japanese technology, Unolok needles are ground on high-precision machines; assembled and packed on totally automated machines. Digital vision inspection cameras ensure 100% automated quality control, enabling superior sharpness to facilitate comfortable insertion.

COMPONENTS

NEEDLE HUBS

Made of non-toxic, medical-grade Polypropylene compatible with the injectable medication color coded to make way for range of sizes.

CANNULA

BURR FREE – 3 FACET LANCET POINT are made from SUS 304 stainless steel

PACKING

The needles are packed in tamper-evident blister packs having easy peelability, medical-grade paper which ensure reliability of seal and product sterility.

PRECAUTIONS

Before use, check the integrity of the sealing of the pack. Store in clean dry & insect free environment. Avoid exposure to excessive heat. Avoid exposure to acidic/corrosive fumes, to avoid corrosion of rust. Dispose in sharps container. Avoid piercing the needle through the wall of hard plastic container, it may result in needle blockage.

COLOUR CODE: As per ISO:6009.



COLOUR CODE AND SPECIFICATION

SIZE		COLOUR	LENGTH							
OD (mm)	Guage		Colour & Code	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
0.40	27G	Grey								
0.45	26G	Brown								
0.50	25G	Orange								
0.55	24G	Medium purple								
0.60	23G	Blue								
0.70	22G	Black								
0.80	21G	Green								
0.90	20G	Yellow								
1.10	19G	Creamish								
1.20	18G	Pink								
1.60	16G	White								

Conform to : ISO: 7864/IS: I O654 standards.

Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.



Scan to
Visit Product Page

