

# KOJAK<sup>®</sup>



## SAFETY BOX

For used syringes, needles & sharps



**MDSAP**  
MEDICAL DEVICE SINGLE AUDIT PROGRAM

# FEATURES/BENEFITS:

- Easy Assembly without training
- Direction of use on the box in multilingual text as well as pictorial directions and sequenced code on flaps.
- Handle in the center avoids spillage/tilted Box during transportation.
- For Disposal -As per country specific waste disposal national policy guidelines.
- Seal edges with Adhesive Surgical Tape in case Box is intended to be used to collect Surgical Blades.
- PQS No. E010/006<sup>\*\*</sup>
- Conforms to WHO standard WHO/PQS/E I O/SBO 1.1<sup>\*\*</sup>

## CARDBOARD SAFETY BOX SPECIFICATION:

Volume	Material	Thickness of walls	Grammage
5 ltr / 10 ltr.	Yellow Duplex Kraft Board	Over 1.10mm	Over 850gsm

## TECHNICAL DIMENSIONS

Weight Fully Assembled	275 gms (min) (5 ltr)	450 gms (min) (10 ltr)
Diameter of syringe insert hole	38 mm (5 ltr)	42mm (10 ltr)
Cartons	5 Ltr. (25 pcs)	10 Ltr. (20 pcs)
Shipping Weight	5 Ltr 7.5 kg (min)	10 Ltr 10 kg (min)

## CAPACITY

S. No.	Size of No. Syringe	With Needle No. Of Syringes Without Agitation		*Without Needle No. Of Syringes Without Agitation		With Needle No. Of Syringes Without Agitation		Without Needle No. Of Syringes Without Agitation	
		5 ltr	10 ltr	5 ltr	10 ltr	5 ltr	10 ltr	5 ltr	10 ltr
	Box Size	5 ltr	10 ltr	5 ltr	10 ltr	5 ltr	10 ltr	5 ltr	10 ltr
1.	0.5ml/2ml	100	200	275	630	150	275	295	650
2.	5ml	65	120	140	320	120	190	170	360

\* Defanged/Cut Needle (5 ltr - suitable for single use syringes 1 ml to 10 ml; 10 ltr suitable for single use syringes 1 ml to 20 ml)