

UNOLOK[®]

STERILE SINGLE USE NEEDLES



INDUSTRIAL SYRINGES
& MEDICAL DEVICES, LTD.
ADDRESS: 100, 10445, 2019

UNOLOK[®]

STERILE SINGLE USE NEEDLES

Manufactured using premium-quality cold-rolled stainless steel strips with the latest Japanese technology, Unolok needles are ground on high-precision machines: assembled and packed on totally automated machines. Vision cameras ensure 100% automated quality control, enabling superior sharpness to minimize pain.

COMPONENTS

NEEDLE HUBS

Made of non-toxic, medical-grade Polypropylene compatible with the injectable medication.

NEEDLE COVER

Made of non-toxic, medical-grade Polypropylene

CANNULA

BURR FREE - 3 FACET LANCET POINT are made from SUS 304 stainless steel

PACKING

The needles are packed in tamper-evident blister packs having easily peelable, medical-grade paper which ensure reliability of seal and product sterility.

PRECAUTIONS

Before use, check the integrity of the sealing of the pack. Store in clean Dry & insect free place. Avoid exposure to excessive heat. Avoid exposure to acidic/corrosive fumes, to avoid corrosion/rust.

Avoid piercing the needle through the wall in hard plastic container, it may result in needle blockage.

COLOUR CODE: As per ISO:6009.



COLOUR CODE AND SPECIFICATION

SIZE		COLOUR	LENGTH							
OD (mm)	Guage		Colour & Code	3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
0.40	27G	Grey								
0.45	26G	Brown								
0.50	25G	Orange								
0.55	24G	Medium purple								
0.60	23G	Blue								
0.70	22G	Black								
0.80	21G	Green								
0.90	20G	Yellow								
1.10	19G	Creamish								
1.20	18G	Pink								
1.60	16G	White								

Conform to : ISO: 7864/IS: 10654 standards.

Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.



Scan to
Visit Product Page

