


**techno[®]
cut**

GLASS VAN[®]

Scalpel & Surgical Blades

The edge of Perfection




INDUSTRIAL DIVISION
A MEDICAL DEVICE LISTED
APRIL 20 1998 2016

- ◆ Ideal Choice for Surgeons, Dentists, Chiropodists
- ◆ Made from premium quality Swiss Carbon Steel & Stainless Steel
- ◆ Manufactured as per ISO : 7740
- ◆ US FDA Regulated
- ◆ Individually packed in an aluminium peel pack which is protected by VCI coated paper Sterility with Gamma Radiation of 2.5 Mrad to ensure sterility assurance level of 10^{-6} whilst preventing corrosion

Product Range

Category	Commercial Designation No.
* For small-size fittings (Handle No.3)	10, 10A, 11, 12, 12D, 13, 15, 15C
* For large-size fittings (Handle No.4)	18, 19, 20, 21, 22, 23, 24, 25
Miniature Blades	6100, 6200, 6400, 6700, Minimal Incision Blade No. 67
Podiatry Blades	312, 313, 314, 316
Stitch Cutters	Stitch Cutter



Packing Dimensions for blades

Size of Shipping Carton (mm)	Primary Packaging (pcs)	Shipping Carton (pcs)	Gross Weight (Kg.)	Volumetric Weight (Kg.)
*540 x 180 x 180	100	5000	10.80	3.00

