


# UNOLOK<sup>®</sup>

## STERILE SINGLE USE INSULIN SYRINGES



Unolok Insulin Syringes are manufactured on totally automated printing, assemble and packaging lines as per the requirements of the most stringent quality standards. All the needles are siliconised to minimise pain.

Presence of low dead space avoids wastage of medication and ensures accurate dosage.



INDIA: PVT. LTD. PUNJAB  
& MEDICAL SERVICES, LTD.  
ADDRESS: 110041 2011

# COMPONENTS

## BARRELS, PLUNGERS, NEEDLE COVERS AND BOTTOM CAPS (MULTIPACKS)

Are made of non-toxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application..

## NEEDLE

It comes integrated with the barrel.

## GASKETS

They are made of natural Rubber/Thermoplastic Elastomer are chemically inert and compatible with all kinds of medications for short term contact drug delivery application

## PRINTING

Graduation scale is available with U 40 and U 100 versions.

1ml is available in both U-40 and U-100, 0.3ml available in U-100

## PACKAGING

Insulin syringes are unit packed in a double laminated plastic film/blister pack.

Insulin syringes are also available in economic multipack of 10 units in a double-laminated plastic film.

Colour Code	 RED	 ORANGE	 RED	 ORANGE	 RED	 ORANGE	 RED	 ORANGE	 RED	 ORANGE	 RED	 ORANGE
-------------	---	--	---	--	---	--	---	--	---	--	---	--

	U-40	U-100	U-40	U-100	U-40	U-100	U-40	U-100	U-40	U-100	U-40	U-100
Needle Size	29G x 1/2" 0.33 x 13mm		29G x 1/2" 0.33 x 13mm		30G x 5/16" 0.30 x 8mm		30G x 5/16" 0.30 x 8mm		*31G x 3/16" 0.25 x 5mm		*31G x 3/16" 0.25 x 5mm	
	0.5/1 ml		0.5/1 ml		0.5/1 ml		0.5/1 ml		1 ml		1 ml	
Description	Multipack		Unitpack		Multipack		Unitpack		Multipack		Unitpack	