


**techno<sup>®</sup>  
cut**

**GLASS VAN<sup>®</sup>**

# Scalpel & Surgical Blades


*The edge of Perfection*



- ◆ Ideal Choice for Surgeons, Dentists, Chiropodists
- ◆ Made from premium quality Swiss Carbon Steel & Stainless Steel
- ◆ Manufactured as per ISO : 7740
- ◆ US FDA Regulated
- ◆ Individually packed in an aluminium peel pack which is protected by VCI coated paper Sterility with Gamma Radiation of 2.5 Mrad to ensure sterility assurance level of 10<sup>-6</sup> whilst preventing corrosion

## Product Range

Category	Commercial Designation No.
* For small-size fittings (Handle No.3)	10, 10A, 11, 12, 12D, 13, 15, 15C
* For large-size fittings (Handle No.4)	18, 19, 20, 21, 22, 23, 24, 25
Miniature Blades	6100, 6200, 6400, 6700, Minimal Incision Blade No. 67
Podiatry Blades	312, 313, 314, 316
Stitch Cutters	Stitch Cutter



## Packing Dimensions for blades

Size of Shipping Carton (mm)	Primary Packaging (pcs)	Shipping Carton (pcs)	Gross Weight (Kg.)	Volumetric Weight (Kg.)
*540 x 180 x 180	100	5000	10.80	3.00



Scan to Visit Product Page