



DISPO VAN[®]

INSULIN SYRINGE

Making Diabetes Management Affordable

Now
introducing
31G with 3/16" (5mm)
extra fine needle



Product Range

Needle Size	Colour code		Description
29G x 1/2"	U-40	RED	Multipack
0.33 x 13mm	U-100	ORANGE	
29G x 1/2"	U-40	RED	Unit Pack
0.33 x 13mm	U-100	ORANGE	
30G x 5/16"	U-40	RED	Multipack
0.30 x 8mm	U-100	ORANGE	
30G x 5/16"	U-40	RED	Unit Pack
0.30 x 8mm	U-100	ORANGE	
31G x 15/64"	U-40	RED	Multipack
0.25 x 6mm	U-100	ORANGE	
31G x 15/64"	U-40	RED	Unit Pack
0.25 x 6mm	U-100	ORANGE	
31G x 3/16"	U-40	RED	Multipack
0.25 x 5mm	U-100	ORANGE	
31G x 3/16"	U-40	RED	Unit Pack
0.25 x 5mm	U-100	ORANGE	

Warning

- Do not share syringes to avoid cross contamination.
- Do not attempt to re-sterilise or reuse.

Precautions

- Before use check the integrity of tamper proof heat-stacked sealing where the end-caps and barrel meet (multi pack). Twist off the bottom cap and the needle cover to expose the plunger and the needle respectively.
- Pull straight off to avoid damaging needle point.
- Manufactured as per ISO:8537



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.

174/25, Ballabgarh, Faridabad (Haryana), Pincode-121004. India | www.hmdhealthcare.com | Info@hmdhealthcare.com |



Scan to
Visit Product Page
Revision No:(Jun. 2024/Rev.-01)