


# techno<sup>®</sup> **PRO** cut



## Cartridge - Blade Management System

*"ALL IN ONE" ATTACHES. DETACHES. PROTECTS.*



- Ideal Choice for Surgeons, Dentists, Chiropractors
- Made from premium quality Swiss Stainless Steel
- Manufactured as per ISO : 7740
- US FDA Regulated
- Allows for effortless, completely safe attachment of standard surgical blade to No3 stainless steel scalpel handle.
- Features a protective sliding sheath ensuring accident free handling.
- Easy to detach blade contained in cartridge for safe disposal.
- Individually packed in an thermoformed blister pack which is protected by VCI coated paper Sterility using ethylene oxide gas to ensure sterility assurance level of 10<sup>-6</sup>.

## Blade Management System

MS Blade Management System “all-in-one” cartridge system provides the weight and feel of a traditional scalpel handle with the cartridge device that uses a standard scalpel handle that attaches directly to the surgical blade.



## Product Range

Category	Commercial Designation No.
* For small-size fittings (Handle No.3)	10, 11, 12, 12D, 15, 15C
* For large-size fittings (Handle No.4)	20, 21, 22, 23, 24, 25

## Packing Dimensions for blades

Size of Shipping Carton (mm)	Primary Packaging (pcs)	Shipping Carton (pcs)	Gross Weight (Kg.)	Volumetric Weight (Kg.)
*555 x 225 x 255	10	1000	5.20	3.00



Scan to Visit Product IFU Page