



UNOLOK[®]

STERILE SINGLE
USE SYRINGES



FEATURES

Unolok Syringes and Needle comprises of a Luer Lock hub, for a more secure and better fitment.

Owing to superior grades of Polypropylene, high clarity plunger and barrel are molded for better visibility of the medicament.

A siliconized rubber piston ensures smooth injection and aspiration motion without scratching the plastic barrel.

The in-house manufactured cannula with three facet lancet pointed tip is well siliconized to reduce friction resulting in smooth penetration.

Unolok Syringes are preferred for usage by dentists and cardiologists as they require a more steady and secure solution to meet their requirements.

PRODUCT RANGE

Capacity	Type of Nozzle	Standard Needle Dimensions
1 ml	Luer mount/luer lock with central nozzle	26G x 1/2" 0.45 x 13
2ml	Luer mount/luer lock with central nozzle	24G x 1" 23G x 1" 23G x 1 1/2" 22G x 1" 0.55 x 25 0.60 x 25 0.60 x 38 0.70 x 25 22G x 1 1/4" 22G x 1 1/2" 21G x 1" 21G x 1 1/2" 0.70 x 32 0.70 x 38 0.80x 25 0.80 x 38
3ml	Luer mount/luer lock with central nozzle	Specifications same as 2ml
5ml	Luer mount/luer lock with central nozzle	Specifications same as 2ml
10ml	Luer mount with eccentric nozzle Luer lock with central nozzle	22G x 1 1/4" 22G x 1 1/2" 21G x 1 1/2" 20G x 1 1/2" 0.70 x 32 0.70 x 38 0.80 x 38 0.90 x 38 Same as above
20m	Luer mount with eccentric nozzle Luer lock with central nozzle	Specifications same as 10ml
50m	Luer mount with eccentric nozzle Luer lock with central nozzle	Without needle

Conform to: ISO 7886-1/1S- 10258 (Part-1) Standards

Also available in B.P

Caution

Reuse of product may lead to cross contamination/ infection.

Store in clean, dry and insect free place

Avoid exposure to acidic/corrosive fumes to avoid corrosion / rust

Warning

Do not use, if the pack is previously open or damaged

Use after checking the Expiry Date

Do not attempt to re-sterilize and re-use

Shelf life

The Shelf Life of the product is 5 years from manufacturing



Scan to Visit Product Page



Revision No: Aug. 2024/Rev.-01

HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.

174/25, Ballabgarh, Faridabad (Haryana), Pincode-121004. India | www.hmdhealthcare.com | Info@hmdhealthcare.com |