




# dispojekt<sup>TM</sup>

SINGLE USE SAFETY (S.I.P.) NEEDLE



Safety from Accidental Needle Injuries  
taking care of those who take  
care of us

---

# Features of dispojekt SAFETY SIP SHIELD

- With safety cover
- Translucent safety mechanism activated by one hand.
- Audible click sound confirms locking of safety device
- Does not require a change in injection technique
- Sharp three –beveled needle to minimize discomfort
- Siliconized needle minimized penetration and gliding force
- Colour coded hub to conform to ISO 6009
- Available in 21G, to 26G in 1 inch and 1.5 inch.
- Needle are offered in thin wall bevel for sizes 21G to 25G, while the 26G option comes with a regular wall bevel design

COLOR CODE		BROWN	ORANGE	MEDIUM PURPLE	BLUE	BLACK	GREEN	
SIZE	OD (MM)	0.45	0.50	0.55	0.60	0.70	0.80	
	GAUGE	26G	25G	24G	23G	22G	21G	
LENGTH		1" ●●●●●			1½" ●●●●●			

## Intended use

Used with or without needle luer mount syringes for single use hypodermic / intravascular / intramuscular injection.

## Packaging details

The Dispojekt single use Safety Needle with SIP (Sharps Injury Prevention) Shield are individually packed in transparent thermoformed peelable blister pack as a strip of five pieces.

Conform to ISO 7864/IS 10654  
 Conform to: ISO 23908  
 (for sharp injury prevention feature)



### Warning

- Do not use if the pack is previously opened or damaged.
- Do not attempt to re-sterilize and reuse.

### Caution

- Reuse of product may lead to cross contamination / infection. Store in clean, dry & insect free place.
- Avoid exposure to acidic / corrosive fumes to avoid corrosion /rust.

### Storage conditions:

- Store in clean, dry & insect free place. Avoid exposure to acidic / corrosive fumes to avoid corrosion /rust.



Scan for Product Film



Scan to watch the tutorial



Scan for 2D Video