



DISPO VAN[®]

STERILE SINGLE USE NEEDLES



STERILE SINGLE USE NEEDLE

- Manufactured using premium-quality cold-rolled stainless steel strips with the latest Japanese Technology.
- DispoVan needles are assembled and packed on totally automated machines. Vision cameras ensure 100% automated quality control.
- Needle Hub made of non-toxic, medical-grade Polypropylene compatible with the injectable medication.
- Cannula Burr Free 3 Facet lancet point with siliconized needle ensures superior sharpness for better patient comfort.

COLOUR CODE AND SPECIFICATION

SIZE		COLOUR		LENGTH						
OD (mm)	Gauge	Colour & Code		3/8"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
0.30	30G	Yellow	Yellow		Yellow					
0.40	27G	Grey	Grey		Grey		Grey			
0.45	26G	Brown	Brown	Brown	Brown			Brown		Brown
0.50	25G	Orange	Orange			Orange		Orange		Orange
0.55	24G	Medium purple	Medium purple					Medium purple	Medium purple	Medium purple
0.60	23G	Blue	Blue					Blue	Blue	Blue
0.70	22G	Black	Black				Black	Black	Black	Black
0.80	21G	Green	Green					Green	Green	Green
0.90	20G	Yellow	Yellow					Yellow	Yellow	Yellow
1.10	19G	Creamish	Creamish						Creamish	Creamish
1.20	18G	Pink	Pink				Pink	Pink		Pink
1.60	16G	White	White							

Conform to : ISO: 7864/IS: 10654 standards.

Before use check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.

174/25, Ballabgarh, Faridabad (Haryana), Pincod-121004. India | www.hmdhealthcare.com | Info@hmdhealthcare.com |



Scan to
Visit Product Page

Revision No: Sep. 2024/Rev-01