



DISPOWAY

**3 WAY
STOP COCK**




Features :

- *Lipid resistant*
- *With two female and one male luer lock port with caps*
- *Leak proof*
- *Secure fitment*
- *Transparent flow channel*
- *Confirms the movement to desired location*

Intended Use :

3 Way Stop Cock are used for intravenous blood transfusion or the transfusion of blood derivatives and infusion of nutritive or medicinal fluids. 3 Way Stop Cock can be used for infusion of more than one fluid at a single time.



Caution

Reuse of product may lead to cross contamination / infection. The use of this product is restricted to qualified medical personnel. Store in clean, dry & insect free place. Avoid exposure to acidic / corrosive fumes.

Contraindication

Product should not be used in patient with known hypersensitivity to any of material used. (Material used : Poly propylene , polycarbonate & Poly ethylene)

Shelf life

The shelf life of the product is 5 years from the date of manufacturing.

