



KOJAK[®] SELINGE

THE KOJAK AUTO-DISABLE HYPODERMIC SYRINGE

To cater to your need for a safe injection



KEY BENEFITS OF KOJAK SELINGE:

Mitigates the risk of cross-contamination and blood borne diseases

The fixed needle ensures that the needles can't be used again.

Minimizes healthcare costs at macro-level

Available in detachable and fixed needle variants

For the curative range of products reconstitution of medications is possible

Available for both curative and fixed dose immunization

Capacity ranges from 0.1 ml to 20ml

For Immunisation Product Range

Size	0.1 ml		0.2 ml	0.25 ml	0.3 ml	0.5 ml
Needle Size	26G x 3/8" 0.45x10	27G x 3/8" 0.40x10	23G x 1" 0.60x25	23G x 1" 0.60x25	23G x 1" 0.60x25	23G x 1" 0.60x25
WHO PQS No.	E008/045		E008/092	E008/095	E008/084	E008/007 Two Ring E008/039 Three Ring

· Conform to: ISO: 7886-3, IS 10258 (Part-3)

· Auto Disable Syringe with Fixed Needle

For Curative Product Range

Size	1 ml		2ml	5ml	10ml	20ml
Needle Size	26G X 1/2" 0.45x13	24G X 1" 0.55x25	26G x 1/2" 0.45x13	23G x 1" 0.60x25	22G x 1 1/4" 0.70x32	21G x 1 1/2" 0.80x38
WHO PQS No.	23G X 1" 0.60x25	22G X 1" 0.70x25	E013/001	E013/018	E013/020	E013/039
	E013/056					

· Conform to: ISO: 7886-4

· Auto Disable Syringe with/without-Fixed/Detachable Needle

Caution

Don't push the plunger before taking the vaccine/medication.

Insert the needle in the vial and pull the plunger slowly until you notice restriction. Do not pull the plunger further after restriction

Shelf life

The Shelf Life of the product is 5 years from manufacturing



Scan to View Product Video



Scan to Visit Product Page

Revision No.(Sep. 2024/Rev-01)



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.

174/25, Ballabgarh, Faridabad (Haryana), Pincode-121004. India | www.hmdhealthcare.com | Info@hmdhealthcare.com |