

# UNOLOK<sup>®</sup> ‡

## SCALP VEIN INFUSION THERAPY EQUIPMENT





# PRODUCT RANGE

Needle Gauge	Colour Code	Needle Length		Needle Diameter (mm)
		(mm)	(Inches)	
18G	Pink	19	3/4"	1.20
19G	Cream	19	3/4"	1.10
20G	Yellow	19	3/4"	0.90
21G	Green	19	3/4"	0.80
22G	Black	19	3/4"	0.70
23G	Light Blue	19	3/4"	0.60
24G	Purple	19	3/4"	0.55
25G	Orange	19	3/4"	0.50
26G	Brown	19	3/4"	0.45
27G	Grey	19	3/4"	0.40

## COMPONENTS

### LUER LOCK CONNECTOR

Luer connector are transparent which allow visual monitoring of the flow of the medicaments. Luer adaptors have a 6% luer taper as per ISO: 80369-7 which provides a leak free mating with an infusion set or a syringe of any make and size to prevent accidental leakage

### WINGS

The soft wings conform well to the skin countours and make immobilisation simple.

### LINK TUBE

The link tube is translucent for easy observation of the contents.

### NEEDLE

The needle is thin walled which allows a larger internal diameter and is ideal for rapid infusion and dosage with high viscosity medications. The needle has a short bevel and three facet point to minimise the chances of cross perforation of vein during insertion. Needles are siliconised to reduce friction & pain (for smooth penetration and withdrawal)

### PROTECTIVE SHEILD

In **UNOLOK**<sup>®</sup> assists in ensuring there is no secondary needle sticks and no reuse of device.

### PACKING

Individually packed in transparent thermoformed peelable blister pack

UNIT PACK	SHELF CARTON	SHIPPING CORRUGATED BOX
Ten blister of <b>UNOLOK</b> <sup>®</sup> Scalp Vein sets form a strip. Each unit is easily separable from the strip along with perforated slitting	100 of units of <b>UNOLOK</b> <sup>®</sup> Scalp Vein sets (16 strip of 6 units + 4 units) are packed in a printed duplex board box	10 shelf cartons are packed in 5 ply corrugated cartons. Pack Size 360x325x490mm Gross Wt: 6.5 Kg (approx)



Scan to View the product Video



Scan to Visit product Page

Revision No: (Jun. 2024/Rev-01)



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.

174/25, Ballabgarh, Faridabad (Haryana), Pincode-121004. India | [www.hmdhealthcare.com](http://www.hmdhealthcare.com) | [Info@hmdhealthcare.com](mailto:Info@hmdhealthcare.com) |