



Cathula[®]†



Safety I.V. Cannula with **tipclip**[™]

Taking care of those who take care of us



As part of your job and your passion, you care for hundreds of patients every day. Doctor, nurses and medical staff are depended upon in time of need. As a health professional you can now depend upon someone as well.

WHENEVER YOU ARE ON DUTY, CATHULA + IS ALSO ON DUTY

A protection device  encapsulates the tip of the used needle when introducer needle is extracted from the Catheter; Post Cannulation.  prevents accidental needle stick injury provided that general effective and safe working procedures and precautions are maintained during its use and disposal

AUTO PROTECTION

The hassel-free and auto protection technology of Cathula+ I.V.Cannula is designed to provide 'Involuntary Activated Safety'-an intuitive mechanism that requires no additional training of the health staff and strives to make I.V. insertion and introducer needle extraction completely safe.

INTRODUCER NEEDLE

Siliconised electrofinished, back cut bevelled stainless steel needle.

Instant Flash Back Notch on 24G and 26G Cannula.



Sharp Injury
Prevention Clip.



NEEDLE HUB

Allows blood flash back visualisation.

Thumb rest for proper grip during insertion.

LUER LOCK CAP

Closure to wing housing permits patient mobility when infusion temporarily discontinued.

BLOOD STOPPER

Designed to expedite blood flash back for enabling correct and efficient cannulation.

BLISTER PACKAGE

Cannula is packed in rigid blister and sealed in peel open type medical grade tyvek.

CATHETER

Available in FEP radio detectable with 4 strips

INJECTION PORT CAP

It facilitates extra medication by preventing its back flow.



COMPONENTS



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
174, 178/25 BALLABGARH FARIDABAD, HARYANA
121004 INDIA

Manufacturing Site :
HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
JHARSANTLY-JAJRU ROAD SECTOR-59 PART-II
BALLABGARH FARIDABAD, HARYANA 121004 INDIA

www.hmdhealthcare.com | info@hmdhealthcare.com

EC REP

Advena Ltd.

Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

CE
0297



Scan to
Visit Product Page

Revision No: (Jan.2025/Rev.-01)