



**Cathy<sup>®</sup>** 

**Safety I.V Cannula**  
with **SIPCLIP<sup>™</sup>**

*Taking care of those who take care of us*



As part of your job and your passion, you care for hundreds of patients every day. Doctors, nurses and medical staff are dependent upon in times of need. As a healthcare professional you can now depend on someone as well.

WHENEVER YOU ARE ON DUTY, CATHY+ IS ALSO ON DUTY

A protection device encapsulates the tip of the used needle when introducer needle is extracted from the catheter Post Cannulation. SIPCLIP prevents accidental needle stick injury provided that general effective and safe working procedures and precautions are maintained during its use and disposal.

AUTO PROTECTION

The hassle-free and auto protection technology of CATHY+ I.V Cannula is designed to provide “Involuntary Activated Safety” – an intuitive mechanism that requires no additional training of the health staff and strives to make I.V insertion and introducer needle extraction completely safe.

COMPONENTS

NEEDLE

Siliconised electrofinished, back cut bevelled stainless steel needle  
Instant Flash Back Notch on 24G and 26G Cannula



LUER LOCK CAP

Closure to wing housing permits patient mobility when infusion temporarily discontinued.

BLOOD STOPPER

The stopper has fine air vent to expedite blood flash back for enabling correct cannulation.

NEEDLE HUB

Allows blood flash back visualisation.  
Thumb rest for proper grip during insertion.

CATHETER

Available in VOLEX™ polyurethane radio detectable.

INJECTION PORT CAP

It facilitates extra medication and also prevents its back flow.

BLISTER PACKAGE

Cannula is packed in rigid blister and sealed in peel open type medical grade paper

PRODUCT RANGE

Colour Code	Gauge	Catheter Ext. Dia x Length(mm)	Water flow-rate (ml/min)
Orange	14G	2.2 x 45	310
Grey	16G	1.7 x 45	200
White	17G	1.5 x 45	140
Green	18G	1.2 x 38	105
		1.2 x 45	100
Pink	20G	1.0 x 32	64
Blue	22G	0.8 x 25	38
Yellow	24G	0.7 x 19	22
Violet	26G	0.6 x 19	15



Scan to  
Visit Product Video



Scan to  
Visit Product Page

