

# dispojekt<sup>TM</sup>

SINGLE USE SAFETY (S.I.P.) NEEDLE



Safety from Accidental Needle Injuries  
taking care of those who take  
care of us

---

# Features of dispojekt SAFETY SIP SHIELD

- With safety cover
- Translucent safety mechanism activated by one hand.
- Audible click sound confirms locking of safety device
- Does not require a change in injection technique
- Sharp three –beveled needle to minimize discomfort
- Siliconized needle minimized penetration and gliding force
- Colour coded hub to conform to ISO 6009
- Available in 21G, to 26G in 1 inch and 1.5 inch.
- Needle are offered in thin wall bevel for sizes 21G to 25G, while the 26G option comes with a regular wall bevel design

## COLOR CODE



SIZE	OD (MM)	0.45	0.50	0.55	0.60	0.70	0.80
	GAUGE	26G	25G	24G	23G	22G	21G
LENGTH		1"      			1½"      		

## Intended use

Used with or without needle luer mount syringes for single use hypodermic / intravascular / intramuscular injection.

## Packaging details

The Dispojekt single use Safety Needle with SIP (Sharps Injury Prevention) Shield are individually packed in transparent thermoformed peelable blister pack as a strip of five pieces.

Conform to ISO 7864/IS 10654

Conform to: ISO 23908

(for sharp injury prevention feature)



### Warning

- Do not use if the pack is previously opened or damaged.
- Do not attempt to re-sterilize and reuse.

### Caution

- Reuse of product may lead to cross contamination / infection. Store in clean, dry & insect free place.
- Avoid exposure to acidic / corrosive fumes to avoid corrosion / rust.

### Storage conditions:

- Store in clean, dry & insect free place. Avoid exposure to acidic / corrosive fumes to avoid corrosion / rust.



Scan for Product Film



Scan to watch the tutorial



Scan for 2D Video



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.  
A28652 ISO 13485 : 2016



**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
174, 178/25 BALLABGARH FARIDABAD, HARYANA  
121004 INDIA

EC REP

**Advena Ltd.**  
Tower Business Centre, 2nd Flr.,  
Tower Street, Swatara, BKR 4013 Malta