



DISPO VAN[®]

INSULIN SYRINGE

Making Diabetes Management Affordable

Now
introducing
**31G with 3/16" (5mm)
extra fine needle**

COMPONENTS

BARRELS, PLUNGERS, NEEDLE COVERS AND BOTTOM CAPS (MULTIPACKS)

Are made of non-toxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application..

NEEDLE

It comes integrated with the barrel.

GASKETS

They are made of natural Rubber/Thermoplastic Elastomer are chemically inert and compatible with various formulation of insulin.

PRINTING

Graduation scale.

1ml is available in both U-40 and U-100, 0.5ml available in U-100

Colour Code  RED	U-40 1 ml	U-40 1 ml	U-40 1 ml	U-40 1 ml	U-40 1 ml	U-40 1 ml
Needle Size	29G x 15/64" (0.33x6mm)	29G x 15/64" (0.33x6mm)	30G x 15/64" (0.30x6mm)	30G x 15/64" (0.30x6mm)	31G x 3/16" (0.25x5mm)	31G x 3/16" (0.25x5mm)
	29G x 1/2" (0.33x13mm)	29G x 1/2" (0.33x13mm)	30G x 5/16" (0.30x8mm)	30G x 5/16" (0.30x8mm)	31G x 15/64" (0.25x6mm)	31G x 15/64" (0.25x6mm)
Description	Multipack	Unit pack	Multipack	Unit pack	Multipack	Unit pack

Colour Code  ORANGE	U-100 0.5/1 ml	U-100 0.5/1 ml	U-100 0.5/1 ml	U-100 0.5/1 ml	U-100 0.5/1 ml	U-100 0.5/1 ml
Needle Size	29G x 15/64" (0.33x6mm)	29G x 15/64" (0.33x6mm)	30G x 15/64" (0.30x6mm)	30G x 15/64" (0.30x6mm)	31G x 3/16" (0.25x5mm)	31G x 3/16" (0.25x5mm)
	29G x 1/2" (0.33x13mm)	29G x 1/2" (0.33x13mm)	30G x 5/16" (0.30x8mm)	30G x 5/16" (0.30x8mm)	31G x 15/64" (0.25x6mm)	31G x 15/64" (0.25x6mm)
Description	Multipack	Unit pack	Multipack	Unit pack	Multipack	Unit pack



Scan to
Visit Product Page



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
174, 178/25 BALLABGARH FARIDABAD, HARYANA
121004 INDIA

EU REP

Advena Ltd.
Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta