



**DISPO VAN<sup>®</sup>**

**SINGLE USE SYRINGES**

**Most Trusted Brand of India**

# THE STERILE HYPODERMIC SYRINGES FOR SINGLE USE ARE PREFERRED FOR THE FOLLOWING REASONS:

- Manufactured on highly automated machines with vision inspection.
- Prevent the spread of blood-borne infections and illnesses.
- Mass-production, allowing greater accessibility and easy availability, meeting any requirements.
- Affordable and widely available.
- Cannula made with premium stainless steel and engineered for a conical and sharp bevel ensuring comfortable insertion.

Capacity	Type of Nozzle	Standard Needle Dimensions
1 ml	Luer mount/luer lock with central nozzle	26G x 1/2" 0.45 X 13
2ml & 2.5mL	Luer mount/luer lock with central nozzle	24G x 1" 23G x 1" 23G x 1/2" 22G x 1" 0.55 X 25 0.6 X 25 0.6 X 38 0.7 X 25 22G X 1/4" 22G X 1/2" 21G X 1" 21G X 1/2" 0.7 X 32 0.7 X 38 0.8 X 25 0.8 X 38
3ml	Luer mount/luer lock with central nozzle	Specifications same as 2ml
5ml	Luer mount/luer lock with central nozzle	Specifications same as 2ml
10ml	Luer mount with eccentric nozzle Luer lock with central nozzle	22G x 1/4" 22G x 1/2" 21G x 1/2" 20G x 1/2" 0.7 X 32 0.7 X 38 0.8 X 38 0.9 X 38 Same as above
20ml	Luer mount with eccentric nozzle Luer lock with central nozzle	Specifications same as 10ml
50ml	Luer mount with eccentric nozzle Luer lock with central nozzle	Without needle

Conform to: ISO 7886-1/IS- 10258 (Part-1) standards

## Caution

Reuse of product may lead to cross contamination/ infection.

Store in clean, dry and insect free place

Avoid exposure to acidic/corrosive fumes to avoid corrosion / rust

## Warning

Do not use, if the pack is previously open or damaged

Use after checking the Expiry Date

Do not attempt to re-sterilize and re-use

## Shelf life

The Shelf Life of the product is 5 years from manufacturing

Conform to : ISO:7886-1/IS: 10258 (Part-1) standards (Sterile hypodermic syringe for single use).  
Before use check the integrity of the sealing of the pack.  
Do not use if previously opened or damage.



**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
174, 178/25 BALLABGARH, FARIDABAD, HARYANA  
121004 INDIA

Manufacturing Site 2 :  
**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
JHARSAINTRY-JAJRU ROAD SECTOR-59 PART-II  
BALLABGARH, FARIDABAD, HARYANA 121004 INDIA

Manufacturing Site 3 :  
**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
99/25 BALLABGARH, FARIDABAD, HARYANA  
121004 INDIA

**Advena Ltd.**  
Tower Business Centre, 2nd Flr.,  
Tower Street, Swatar, BKR 4013 Malta



Scan to  
Visit Product Page