



3 WAY STOP COCK

With extension tube




Features :

- Three Way Stop Cock with Extension Tube integrated with Kink resistant PVC tubing
- Facilitates drug administration for additional lines together
- With two female and one male luer lock port with caps for secure fitment
- Leak proof body withstands pressure up to 8 bars (115 PSI) when applied in normal conditions
- Uniformity leads to smooth infusion flow with minimal connection.
- Transparent flow channel

Intended Use :

3 Way Stop Cock Extension Tube are used for intravenous transfusion of blood derivatives and infusion of nutritive or medicinal fluids. 3 Way Stop Cock Extension Tube can be used for infusion of more than one fluid at a single time.



Caution

Reuse of product may lead to cross contamination / infection.
The use of this product is restricted to qualified medical personnel.
Store in clean, dry & insect free place.
Avoid exposure to acidic / corrosive fumes.

Shelf life

The shelf life of the product is 5 year from the date of manufacturing.

Contraindication

Product should not be used in patient with known hypersensitivity to any of material used.
(Material used : Polypropylene , polycarbonate & Polyethylene)



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.

99/25, Ballabgarh, Faridabad (Haryana), Pincode-121004. India | www.hmdhealthcare.com | Info@hmdhealthcare.com |

www.hmdhealthcare.com | info@hmdhealthcare.com



Scan for video

Revision No: (Jun. 2025/Rev.-02)