



For used syringes, needles & sharps



FEATURES/BENEFITS:

- Easy Assembly without training
- Direction of use on the box in multilingual text as well as pictorial directions and sequenced code on flaps.
- Handle in the center avoids spillage/tilted Box during transportation.
- For Disposal -As per country specific waste disposal national policy guidelines.
- Seal edges with Adhesive Surgical Tape in case Box is intended to be used to collect Surgical Blades.
- WHO PQS Code. E010/006^{☆☆}
- Conforms to WHO standard WHO/PQS/E10/SB01.1^{☆☆}

CARDBOARD SAFETY BOX SPECIFICATION:

Volume	Material	Thickness of walls	Grammage
5 ltr	Yellow Duplex Kraft Board	Over 1.0mm	Over 850gsm
10 ltr	White Duplex Kraft Board	Over 1.0mm	Over 850gsm

TECHNICAL DIMENSIONS

Weight Fully Assembled	275 gms (min) (5 ltr)	450 gms (min) (10 ltr)
Diameter of syringe insert hole	38 mm (5 ltr)	42mm (10 ltr)
Cartons	5 Ltr. (25 pcs)	10 Ltr. (20 pcs)
Shipping Weight	5 Ltr 7.5 kg (min)	10 Ltr 10 kg (min)

CAPACITY

S. No.	Size of No. Syringe	With Needle No. Of Syringes Without Agitation		*Without Needle No. Of Syringes Without Agitation		With Needle No. Of Syringes With Agitation		Without Needle No. Of Syringes With Agitation	
		5 ltr	10 ltr	5 ltr	10 ltr	5 ltr	10 ltr	5 ltr	10 ltr
1.	0.5ml/2ml	100	200	275	630	150	275	295	650
2.	5ml	65	120	140	320	120	190	170	360

* Defanged/Cut Needle (5 ltr - suitable for single use syringes 1 ml to 10 ml; 10 ltr suitable for single use syringes 1 ml to 20 ml)

☆☆ 5 liter safety box



HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
174/25, Ballabgarh, Faridabad (Haryana), Pincode-121004. India | www.hmdhealthcare.com | Info@hmdhealthcare.com |



UL
UL LISTED
ELECTRIC & ELECTRONIC EQUIPMENT
EQUIPMENT