

**DISPO VAN<sup>®</sup>**



# STERILE SINGLE USE ORAL TIP SYRINGES



FOR ORAL / ENTERAL USE ONLY

## FEATURES

- The transparent barrel, crafted from non-toxic medical-grade polypropylene, enhances visibility, and ensures compatibility with various medications.
- The centric oral cone in oral tip syringes offers precise medication delivery and enhances user control and comfort during administration
- The non-luer tip design eliminates the 6% luer taper, preventing potential use with hypodermic needle.
- Sterilized by the ethylene oxide gas to achieve Sterility Assurance Level  $10^{-6}$ .
- The interior surface of the syringe is lubricated with pharmaceutical grade silicon (Polydimethylsiloxane) for easy and smooth piston movement ensures accurate dosing.

## PRODUCT RANGE

Capacity	Smallest Graduation	Minimum Numbered Graduation	Graduated Scale
1 ml	0.25 ml	0.25 ml	1 ml
3 ml	0.1 ml	0.5 ml	3 ml
6 ml	0.2 ml	1 ml	6 ml
10 ml	0.25 ml	1 ml	10 ml

## PACKING AND SHIPPING

Syringe Capacity	Size (mm) (LXWXH)	Gross weight (kg) approx	Quantity per Shipping Pack
1 ml (blister pack)	440x320x470	9.2	20 shelf pack in a carton
3 ml (blister pack)	505x400x310	10.0	12 shelf pack in a carton
6 ml (blister pack)	810x415x490	21.5	24 shelf pack in a carton
10 ml (blister pack)	475x400x440	11.6	18 shelf pack in a carton

- Standards:** Company specification based upon ISO 7886-1
- Shelf life:** The shelf life of product is 5 year from the date of manufacturing.
- Precautions:** Before use check the integrity of the sealing of the pack.
- Warning:** Do not use if the pack is previously opened or damaged.  
Do not attempt to re-sterilize and reuse.  
Discard in a suitable container.

