



PERFEKTUM.TPE

Sterile Single Use Hypodermic Syringe



Trust | Profitable | Economical



PERFEKTUM.TPE

Sterile Single Use Hypodermic Syringe

Features

- PERFEKTUM. TPE brand, is manufactured using fully automated processes. This 3-part syringe, with Gasket, Plunger and Barrel, meets daily physician needs.
- Our manufacturing lines employ high-standard Vision Inspection systems to verify product compliance, checking for Needle Tip, Glue Bead Formation, Needle Bend, and Needle Cover presence in each piece.
- The syringe barrel hold medication or fluid for injection, marked with graduation scaled for accurate dose measurement.
- The tip of the needle has a thin wall and three facet lancet point. Needle is well siliconised to reduce friction for comfortable insertion.

Capacity	Type of Nozzle	Standard Needle Dimensions			
2ml 2.5ml	Luer mount with central nozzle.	24G x 1" 0.55 x 25	23G x 1" 0.60 x 25	23G x 1½" 0.60 x 38	22G x 1" 0.70 x 25
3ml	Luer mount with central nozzle.	22G x 1¼" 0.7 x 32	22G x 1½" 0.70 x 38	21G x 1" 0.8 x 25	21G x 1½" 0.8 x 38
5ml	Luer mount with central nozzle.				
10ml	Luer mount with eccentric nozzle.	24G x 1" 0.55 x 25	23G x 1" 0.60 x 25	23G x 1½" 0.60 x 38	22G x 1" 0.70 x 25
20ml	Luer mount with eccentric nozzle.	22G x 1¼" 0.7 x 32	22G x 1½" 0.70 x 38	21G x 1½" 0.8 x 38	20G x 1½" 0.9 x 38
50ml	Luer mount with eccentric nozzle.	Without Needle			

Conform to: ISO 7886-1/IS: 10258 standards

