

**techno[®]
cut**

GLASS VAN[®]

Scalpel & Surgical Blades

The edge of Perfection



- ◆ Ideal Choice for Surgeons, Dentists, Chiropodists
- ◆ Made from premium quality Swiss Carbon Steel & Stainless Steel
- ◆ Manufactured as per ISO : 7740
- ◆ US FDA Regulated
- ◆ Individually packed in an aluminium peel pack which is protected by VCI coated paper to prevent corrosion sterilised with Gamma Radiation of 2.5 Mrad to ensure sterility assurance level of 10⁻⁶.

Product Range of Blades

Category	Commercial Designation No.
For small-size fittings (Handle No.3)	10, 10A, 11, 12, 12D, 13, 15, 15C, 15CA
For large-size fittings (Handle No.4)	18, 19, 20, 21, 22, 23, 24, 25
Miniature Blades	6100, 6200, 6400, 6700, Minimal Incision Blade No. 67
Podiatry Blades	312, 313, 314, 316
Stitch Cutters	Stitch Cutter



Product Range of Scalpel

Category	For small size fittings (Handle no. 3)	For large size fittings (Handle no. 4)
Commercial Designation No.	10, 11, 12, 12D, 15, 15C	20, 21, 22, 23, 24, 25



NIRAJ INDUSTRIES PVT. LTD.

923, SECTOR-68, IMT Faridabad - 121004, Haryana, India | www.hmdhealthcare.com



Scan to
Visit Product Page

Revision No: Jun. 2025/Rev.-02