

UNOLOK[®] ‡

SCALP VEIN INFUSION THERAPY EQUIPMENT



PRODUCT RANGE

Needle Gauge	Colour Code	Needle Length		Needle Diameter (mm)
		(mm)	(Inches)	
18G	Pink	19	3/4 "	1.20
19G	Cream	19	3/4"	1.10
20G	Yellow	19	3/4"	0.90
21G	Green	19	3/4"	0.80
22G	Black	19	3/4"	0.70
23G	Light Blue	19	3/4"	0.60
24G	Purple	19	3/4"	0.55
25G	Orange	19	3/4"	0.50
26G	Brown	19	3/4"	0.45
27G	Grey	19	3/4"	0.40

COMPONENTS

LUER LOCK CONNECTOR

Luer connector are transparent which allow visual monitoring of the flow of the medicaments. Luer adaptors have a 6% luer taper as per ISO: 80369-7 which provides a leak free mating with an infusion set or a syringe of any make and size to prevent accidental leakage

WINGS

The soft wings conform well to the skin countours and make immobilisation simple.

LINK TUBE

The link tube is translucent for easy observation of the contents.

NEEDLE

The needle is thin walled which allows a larger internal diameter and is ideal for rapid infusion and dosage with high viscosity medications. The needle has a short bevel and three facet point to minimise the chances of cross perforation of vein during insertion. Needles are siliconised to reduce friction & pain (for smooth penetration and withdrawal)

PROTECTIVE SHEILD

In **UNOLOK[®]** assists in ensuring there is no secondary needle sticks and no reuse of device.

PACKING

Individually packed in transparent thermoformed peelable blister pack

UNIT PACK	SHELF BOX	SHIPPING CORRUGATED CARTON
Ten blister of UNOLOK[®] Scalp Vein sets form a strip. Each unit is easily separable from the strip along with perforated slitting	100 of units of UNOLOK[®] Scalp Vein sets (16 strip of 6 units + 4 units) are packed in a printed duplex board box	10 shelf cartons are packed in 5 ply corrugated cartons. Pack Size 360x325x490mm Gross Wt: 6.5 Kg (approx)



Scan to
View the product Video



Scan to
Visit product Page

Revision No:(Jun. 2025/Rev-02)



Manufacturing Site:
HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.
99/25 BALLABGARH FARIDABAD, HARYANA
121004 INDIA