

# UNOLOK<sup>®</sup>

## STERILE SINGLE USE SYRINGES



# FEATURES

Unolok Syringes and Needle comprises of a Luer Lock hub, for a more secure and better fitment.

Owing to superior grades of Polypropylene, high clarity plunger and barrel are molded for better visibility of the medicament.

A siliconized rubber piston ensures smooth injection and aspiration motion without scratching the plastic barrel.

The in-house manufactured cannula with three facet lancet pointed tip is well siliconized to reduce friction resulting in smooth penetration.

Unolok Syringes are preferred for usage by dentists and cardiologists as they require a more steady and secure solution to meet their requirements.

# PRODUCT RANGE


Capacity	Type of Nozzle	Standard Needle Dimensions			
1 ml	Luer mount/luer lock with central nozzle	26G x 1/2" 0.45 x 13			
2ml	Luer mount/luer lock with central nozzle	24G x 1" 0.55 x 25	23G x 1" 0.6 x 25	23G x 1 1/2" 0.6 x 38	22G x 1" 0.7 x 25
		22G x 1 1/4" 0.7 x 32	22G x 1 1/2" 0.7 x 38	21G x 1" 0.8 x 25	21G x 1 1/2" 0.8 x 38
3ml	Luer mount/luer lock with central nozzle	Specifications same as 2ml			
5ml	Luer mount/luer lock with central nozzle	Specifications same as 2ml			
10ml	Luer mount with eccentric nozzle	22G x 1 1/4" 0.7 x 32	22G x 1 1/2" 0.7 x 38	21G x 1 1/2" 0.8 x 38	20G x 1 1/2" 0.9 x 38
	Luer lock with central nozzle	Same as above			
20m	Luer mount with eccentric nozzle Luer lock with central nozzle	Specifications same as 10ml			
50m	Luer mount with eccentric nozzle Luer lock with central nozzle	Without needle			

Conform to: ISO 7886-1/1S- 10258 (Part-1) Standards  
Also available in B.P

**Caution**  
Reuse of product may lead to cross contamination/ infection.  
  
Store in clean, dry and insect free place  
  
Avoid exposure to acidic/corrosive fumes to avoid corrosion / rust

**Warning**  
Do not use, if the pack is previously open or damaged  
  
Use after checking the Expiry Date  
  
Do not attempt to re-sterilize and re-use

**Shelf life**  
The Shelf Life of the product is 5 years from manufacturing

  
**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
174, 178/25 BALLABGARH, FARIDABAD, HARYANA  
121004 INDIA

**Manufacturing Site 2 :**  
**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
JHARSAINLY-JAJRU ROAD SECTOR-59 PART-II  
BALLABGARH, FARIDABAD, HARYANA 121004 INDIA

**Manufacturing Site 3 :**  
**HINDUSTAN SYRINGES & MEDICAL DEVICES LTD.**  
99/25 BALLABGARH, FARIDABAD, HARYANA  
121004 INDIA



Scan to  
Visit Product Page

**EU REP**  
**Advena Ltd.**  
Tower Business Centre, 2nd Flr.,  
Tower Street, Swatar, BKR 4013 Malta